

DNA Sequencing Order Form



Our Address: Sequencing Center, Vaccination Dep., Pasteur Institute of Iran , No. 358, Rajabi Alley, 12 Farvardin St., Tehran, Iran.
E-Mail: pasteurseq@yahoo.com **Fax:** 66496582

Customer Profile:

Name:..... Organization:..... Phone:.....
 Address:
 Mobile:..... E-Mail:..... Fax:.....

Page ...of

Dear Customer, In order to receive successful sequencing results, please fill in this order carefully and send it along with your DNA samples(in 0.5ml tube) to the address above .

Order No.

NO.	Template or Sample Name	Primer Name	Template Type	Template Size Including Insert	DNA Conc. by Spectrophotometer	Template GC% Content	Primer Conc. (pMol)	Primer Tm	Purification* Is needed? Yes/No?
1									
2									
3									
4									
5									
6									
7									
8									
9									
10									
11									
12									
13									
14									

*PCR product needs to be purified. The PCR product may also be purified by us (extra service).

We need the following reagents per reaction: **Primer:** 10µl at 10pMol. **PCR products:** 10µl at 10-100ng/ µl (purified) and 100µl (nonpurified). **Plasmid:** 10µl at 150-300ng/µl. **Cosmid, BAC:** 10µl at 0.5-1.0µg/µl. **Bacterial genomic DNA:** 10 µl at 2-3 µg/µl.